



BLIXTLÁS



2-1 FUNKTION

INSTRUCTIONS FOR USE SOMNA® BLANKET BALANCE

Settle down with the Somna® Blanket Balance. The chains inside the length-wise channels shape the blanket tightly to the body's contours, giving the user a pleasant feeling of being embraced.



USE

In the Somna® Blanket Balance the chains are more centrally placed, in line with the body's position, providing an extra embracing effect.

The Somna® Blanket Balance has a two-in-one function: one chain side and one padded side. Which side is against the body is down to personal choice. Having the chains against the body provides greater tactile stimulation. Try out both sides to see what's best for you. The chain side has a half-moon illustration/icon to easily distinguish it from the padded side.

Use the Somna® Blanket Balance in a duvet cover. You can use the blanket for a whole night's sleep, but we recommend using it for short intervals to begin with, and then increasing gradually. The blanket is so supple that you can also fold it if you require more weight on a certain area of the body. The Blanket is designed to be as cool as possible.

Draping the blanket over the shoulders and letting it adapt to your body shape while seated is nicely relaxing when resting during the daytime.

WASHING INSTRUCTIONS



Maximum 60°C – gentle wash cycle.



No liquids stronger than PCE when dry cleaning



Do not tumble dry



Do not iron



Do not bleach



A zipped blanket should be opened before washing. A laundry bag can be used if desired.

PRODUCT INFORMATION



The Somna® Blanket Balance fulfils SS-EN ISO 12952-2 standards (Ignitability of bedding items – Ignition source: match flame equivalent). Outer fabric in flame-resistant Trevira CS material.

Textile: 100% Polyester Trevira CS, **Filling:** Steel chains & Polyester



The Somna® Blanket Balance is CE approved to the Medical Devices Directive 93/42/EEC.



The Somna® Blanket Balance is a Swedish patented product.

SIZES/WEIGHTS

200 x 140 cm / 8 kg

NOTE

- If the user has reduced lung capacity, seek medical advice before using a weighted blanket.
- Consult a doctor if the user has any surgically inserted devices that may be affected by the steel in the chains.
- Beware the risk of tripping if the blanket is left on the floor.
- Take care if the user has low muscle capacity.
- When using hygiene cover ensure supervision is provided at all times as the cover could result in potential choking hazard.