


Size chart for Somna® Vest

Child

Size	Bust/chest width Min	Max	Waist width Max	Length Centre back
110/122	56	62	68	44
128/140	66	74	80	48
146/158	76	86	92	52



Adult

Size	Bust/chest width Min	Max	Waist width Max	Length Centre back	Ladies' size	Men's size
S	88	98	106	59	38-42 (M-L)	44-48 (S-M))
M	96	110	118	63	44-48 (L-XL)	50-54 (M-L)
L	108	120	128	66	48-52 (XL-XXL)	56-60(XL-XXL)
XL	118	134	142	70	52-56(XXL-XXXL)	60-66(XXL-XXXL)

How to measure the body:

- The chest/bust measurement is taken around the body at chest height, without loosening or tightening the tape measure.
- The waist measurement is taken around the body at navel height, without loosening or tightening the tape measure.
- The vest's waist measurement is adjusted with a belt and can be about 8 cm larger than the bust/chest measurement. The child vest is about 6 cm larger at the waist.

How to measure the garment:

- The chest/bust measurement is taken straight across the garment at the bottom of the armpit when the vest is buttoned up.
- The centre back length is measured from the edge of the neck at the centre of the back, to the bottom edge of the vest.



Measuring the smallest measurement for the size. No seam overlap.



Measuring for the largest measurement for the size. None of the Velcro strip is visible.



Measurement for the length at centre back.

SOMNA AB

Panelvägen 16, 443 61 Stenkullen, Sweden
+46 31 301 17 20
info@somna.se | somna.se

